

# **Introduction to genome tiling microarray analysis**

# Biological motivations

- There are many types of “events” happen at different locations on the genome. For example, protein bindings, epigenetic modifications (DNA methylation and histone modifications), copy number variations, etc.
- It is often of great interests to detect the genomic locations where a specific event happens, or quantify the events along the genome.
- The locations of these events provide explanations for many biological processes.

# An example: transcription factor(TF) binding

- Transcription factors (TF): proteins that binds to specific DNA sequences and control the transcription from DNA to mRNA.
- There are many different types of TFs, each recognize different DNA sequences (motifs).
- The functions of the TFs are important for understanding gene regulatory mechanisms.
- The first step toward the understanding is to detect the TF binding sites (TFBS).

# Traditional Method for Understanding Transcription Regulation




**Gene expression  
microarray analysis**



**Clustering genes by  
expression profile**



**Search conserved sequence  
motifs in cluster promoters**



Very challenging for mammalian genomes!

## Use tiling array (ChIP-chip) to detect TFBS


- Detect genome-wide *in vivo* location of TF and other DNA-binding proteins.
- Can learn the regulatory mechanism of a transcription factor or DNA-binding protein much better and faster.

## Another example: DNA methylation and histone modification

- DNA methylation and histone modifications are chemical modification of DNA molecule.
- The strengths of such modifications varies along the genome, and they are related to gene regulatory and many diseases.
- The methylation or modification strengths can be measured using ChIP-chip or MeDIP-chip.

# Tiling arrays

- The goal is to quantify the events of interests along the genomes, and/or detect the genomic coordinates for the events.
- Work the same as gene expression array (hybridization based), except that the probes are designed to tile up the genome at non-repeat regions.
- Data for probes in the location of interest often behave differently from backgrounds (e.g., bigger intensities).



# Types of tiling arrays

- ChIP-chip: Chromatin ImmunoPrecipitation (ChIP) + tiling array (chip) for detecting transcription factor binding sites or measuring histone modification levels.
- MeDIP-chip: Methyl-DNA ImmunoPrecipitation (MeDIP) + tiling array (chip) for measuring DNA methylation level.
- ArrayCGH (Comparative Genomic Hybridization) for detecting copy number variations.


There's no major differences in array designs. Difference are the ways to prepare biological samples.




# Available platforms for ChIP-chip

	<b># Arrays human genome</b>	<b># Probes / Array</b>	<b># Total Probes</b>	<b>Probe Length</b>	<b>Probe Resolution</b>
<b>Affymetrix</b>	<b>7</b>	<b>6M</b>	<b>42.0M</b>	<b>25mer</b>	<b>35 bp</b>
<b>Nimblegen</b>	<b>10</b>	<b>2.1M</b>	<b>21M</b>	<b>50mer</b>	<b>100 bp</b>
<b>Agilent</b>	<b>21</b>	<b>244K</b>	<b>5.1M</b>	<b>60mer</b>	<b>300 bp in genes; 500 bp in intergenic</b>


# ChIP-chip procedures




# Chromatin ImmunoPrecipitation (ChIP)




# TF/DNA Crosslinking *in vivo*




# Sonication (~500bp)




# TF-specific Antibody



# Immunoprecipitation (IP)




# Reverse Crosslink and DNA Purification






# Amplification




# ChIP-chip Hybridization




# Data from ChIP-chip

- Can be thought as a file with millions of rows and three columns.
  - Each row is for a probe.
  - Columns are chromosome number, probe location on the genome, and signal (intensity values or log fold change).
- To visualize: plot the probe signals against probe locations.



# Identify ChIP-enriched Region

- Controls: sonicated genomic input DNA (non-treated).
- Often 3 ChIP, 3 Ctrl replicates are needed



# ChIP-chip data analysis

- Goal: detect locations of interests (e.g., binding sites, also called “peaks”) based on probe locations and signals.
- Normalization: remove technical artifacts.
- Detection for regions of interests:
  - Many different methods. Fundamentally data from neighboring probes need to be combined to make inference, because the regions of interests often overlap many probes.
  - Easiest method: moving average, then use an arbitrary cutoff.

# Mann-Whitney U-test

- Affy TAS, Cawley et al (*Cell* 2004):
  - Each probe: rank probes signals within [-500bp, +500bp] window.
  - Check whether sum of ChIP ranks is much smaller

	ctrl 1	ctrl 2	ChIP 1	ChIP 2		ctrl 1	ctrl 2	ChIP 1	ChIP 2
probe 1	1.71	2.23	3.02	2.25	probe 1	17	15	13	14
probe 2	4.27	3.10	3.86	4.70	probe 2	6	12	10	3
probe 3	4.06	3.67	4.03	4.74	probe 3	7	11	8	2
probe 4	1.20	0.40	1.31	1.85	probe 4	19	20	18	16
probe 5	4.29	3.95	4.56	4.76	probe 5	5	9	4	1

# TileMap (Ji and Wong, Bioinformatics 2005)

## **STEP 1:**


**Compute a test statistic for each probe to  
summarize probe level information**



## **STEP 2:**


**Combine probe level test statistics of  
neighboring probes to help infer binding regions**

# Probe level test statistic: empirical Bayes approach





# Combining neighboring probes



## TileMap (MA)

1. Compute the probe level test statistic  $t$  for each probe;
2. Compute a moving average statistic to measure enrichment;
3. Estimate FDR.

## TileMap (HMM)

1. Compute the probe level test statistic  $t$  for each probe;
2. Estimate the distribution of  $t$  under  $H_0$  and  $H_1$ ;
3. Model  $t$  by a Hidden Markov Model, and decode the HMM.




# TileMap summary

- Now a part of a software suite CisGenome.
- Windows based GUI.
- Command line version available for Mac and Linux.
- Freely available from:  
<http://www.biostat.jhsph.edu/~hji/cisgenome/>


# CisGenome

(Ji H. et al. Nature Biotechnology, 2008)

## Graphic User Interface




## CisGenome Browser



# MAT: Model-based Analysis of Tiling arrays (Johnson W.E. et al. *PNAS*, 2006)

- Estimate probe behavior by checking other probes with similar sequence on the same array
- Probe sequence plays a big role in signal value.
- Most of the probes in ChIP-chip measures non-specific hybridization.



# Probe Behavior Model




$$\text{Log}(PM_i) = \alpha n_{iT} + \sum_{j=1}^{25} \sum_{k=A,C,G} \beta_{jk} I_{ijk} + \sum_{l=A,C,G,T} \gamma_l n_{il}^2 + \delta \text{Log}(c_i) + \varepsilon_i$$

Baseline on  
number of Ts

A,C,G at each position  
of the 25mer

A,C,G,T Count Square


25mer Copy Number  
along the Genome





By Xiaole Shirley Liu at Harvard

# Probe Standardization

- Fit the probe model array by array
- Divide array probes to bins (3k probes/bin)
- Background-subtraction and standardization (normalization) on a single array;



- Probe signals before and after standardization



# Binding region detection

- Window-based MATscore
  - ChIP without Ctrl

$$MAT(region) = TM(t's \text{ in region}) \sqrt{n_{probe}}$$

- TM: trimmed mean
- Multiple ChIP with multiple Ctrl

$$MAT(region) = \frac{TM(t's \text{ in ChIP}) - TM(t's \text{ in Input})}{\sigma_{Input}} \sqrt{n_{probe}}$$

- More probes, higher t values in ChIP, less variance (fluctuation) → more confident



# To use MAT

- Create a text configuration file (config.txt):

```
[data]
BpmapFolder = /home/bst/student/hwu/Project/Ji/MVHMM/DREAM/rawdata/
CelFolder = /home/bst/student/hwu/Project/Ji/MVHMM/DREAM/rawdata
GenomeGrp =Hs
Group = 111000
[bpmap]
1=Hs_PromPR_v02-3_NCBIV36.bpmap
[cel]
1=IP1.CEL IP2.CEL IP3.CEL CT1.CEL CT2.CEL CT3.CEL
[intensity analysis]
BandWidth =      300
MaxGap =      300
MinProbe =      10
[interval analysis]
Pvalue = 1e-3
```

- Then run “MAT config.txt” at command line.

# MAT summary

- Open source, written in python at <http://chip.dfci.harvard.edu/~wli/MAT/>
- Installation could be tricky.
- Good computational performance.
- Can work with single ChIP, multiple ChIP, and multiple ChIP with controls with increasing accuracy.

# Bioconductor packages for analyzing ChIP-chip data

- Most of the ChIP-chip analysis are done using MAT or CisGenome, so there are relatively fewer R packages.
- Useful ones:
  - rMAT: R implementation of MAT model. Works for Affy ChIP-chip.
  - Ringo (R Investigation of NimbleGen Oligoarrays): works for NimbleGen two-color tiling arrays.
  - Starr: an extension of Ringo, works for Affymetrix arrays.
  - ChIPpeakAnno: annotation of peaks, e.g., find closeby genes, GO terms, DNA sequences, etc.

# rMAT

- R implementation of MAT.
- Works for Affymetrix arrays only.
- Needs bmap file (factory provided file to probe annotations), and raw data file in CEL format.


```
library(rMAT)
bmapFile="Hs.bmap"
### read in data
arrayFile=c("IP.CEL", "Control.CEL")
rawdata=BMAPCelParser(bmapFile, arrayFile, groupName="Sc")

## normalization - run MAT model
normdata=NormalizeProbes(rawdata,method="MAT")

## compute MAT scores and find peaks
RD=computeMATScore(normdata,cName="Control", dMax=600)
Enrich=callEnrichedRegions(RD, dMax=600, dMerge=300,
    nProbesMin=8, method="score", threshold=2)
```

# Copy number variation arrays


- CNV: phenomenon that sections of DNA have abnormal number of copies (deviate from 2).
- Can be detected by SNP arrays (for one sample) or arrayCGH (comparing two samples case vs. control).



# Data from CNV arrays

- Data format are the same, e.g., probe locations and signal, but characteristics are different:
  - No peaks, but long, flat-topped “plateaus”.
  - Heights of the plateaus are discrete, corresponding to different number of copies (integers: 1, 2, 3, ...)

# Example data from arrayCGH





# Analysis of CNV arrays

- Methods are different from ChIP-chip, but still smoothing based to combine neighboring probe information, for example, Hidden Markov Model.

# A list of CNV array software

- Affymetrix:
  - APT: uses a hidden Markov model
  - R package VanillaICE: HMM base. R. Scharpf *et al.* (2008) AOAS
  - R package DNACopy: Circular Binary Segmentation. Olshen *et al.* (2004) Biostatistics
- Illumina:
  - QuantiSNP: S. Colella *et al.* (2007), NAR
  - PennCNV: K. Wang *et al.* (2008), NAR

# Review

- Tiling arrays are DNA microarrays for detecting locational modifications of genome.
- Probes tile up a part of whole genome.
- Still hybridization based (DNA segments stick to probes), same as gene expression arrays.
- Data need to be visualized along genome.
- Location of interests shows some patterns: peaks for TFBS, or plateau for CNV.
- Need to combine data from neighboring probes to make calls.